

## Streamline Patient Care

### Patient Nurse Call Alerts Go Directly to Their Caregiver

Adding the Noti-Fi Smartphone alerting option to the Provider<sup>®</sup> 790 Nurse Call System provides caregivers mobility while still being connected to their patients. When a patient places a call on the system, a text alert instantly routes to their caregiver or caregiver team. A tone or vibration notifies caregivers of an active call while the display shows the room number, type of call priority, and how long the call has been active.


Noti-Fi gives you flexibility and value by using the technology you already have in place. No cell phone data is required because patient call notifications are sent over your existing Wi-Fi LAN. Any Android<sup>\*</sup> phone or tablet can be used for Provider 790 nurse call alerts; facilities can use existing phones/tablets or staff can use their own devices.



# Noti-Fi<sup>™</sup>

### Focus on Patient Safety & Satisfaction by:

- Wireless notifications to caregivers speed response times
- Supports a quiet environment reducing disruptive overhead paging
- Noti-Fi text alerts automatically prioritize by type of call so the most urgent is always displayed first
- Supports transparency and accountability - all staff to patient interactions are logged in the optional reporting package to confirm a timely response
- Overtime reminders automatically alert caregivers when calls have been waiting too long



*\* Android version 6.0 and later*

## With Noti-Fi's Flexible Call Routing, No Call is Missed

The flexibility of Noti-Fi supports **primary caregiver alerting**, **team caregiver alerting**, or any combination of the two.


### Primary Caregiver Alert



Caregivers can back up each other when they respond to patient calls. If the first caregiver is busy and cannot respond, the next caregiver in the team is automatically alerted.

### Team Caregiver Alert

Patient calls simultaneously alert multiple caregiver's smart devices. The closest staff available "Accepts" the call and responds. Other caregivers are notified who has accepted the call. In an emergency situation, all team members can respond.



The intuitive Noti-Fi interface has a full view of active calls, including how long each call has been active

- A tone and/or vibration alerts staff to a new patient call
- Patient calls continue to display on a locked phone or tablet
- Staff manage their On/Off Duty and On Break status directly from their phone or tablet
- "Accepting" a call alerts other caregivers on the team, minimizing steps and stopping multiple people from responding to the call

Jeron's Provider 790 Nurse Call System and Noti-Fi application keeps caregivers mobile and ready to respond to their patients' needs.

Learn more about the Noti-Fi application for Provider 790 Nurse Call, call 800.621.1903 or visit us: [www.jeron.com](http://www.jeron.com)