
Jeron Provider[®] 680 Nurse Call System/ PCS Integration Features

10/01/2004

Features

The following features are provided by Jeron Provider 680 when integrated with the SpectraLink system:

- **Set Nurse Need**
A request for nurse service may be sent from the PCS phone.
- **Set Aide Need**
A request for aide service may be sent from the PCS phone.
- **Send Manual Radio Pocket Page with Optional Tag Message**
- **Call Hunting**
This feature allows a set of handsets to ring in order for a single call. This allows a call that is unanswered at one handset to be forwarded the next handset in the hunt group until the call is answered.
- **Code Groups**
All phones belonging to the code group will ring simultaneously when a code call is placed on the Provider 680 system. The phone will serve as a paging unit only, displaying the call information on the PCS phone.
- **Originate to 680 Station**
Calls can be originated to Provider 680 stations from the PCS phone.
- **Single Touch Answering of Calls**
Calls may be answered by pressing the Start button on the PCS phone.
- **Programmable Number of Rings**
The number of rings can be programmed for the PCS phone.
- **Call Abandonment**
If the PCS phone user fails to terminate the audio connection at the PCS phone, the Provider 680 system will terminate audio on the Provider system for that call when the programmed call abandonment time limit is reached.
- **Priority Ringing**
If a call is placed on the Provider 680 system while the PCS phone is connected in a PBX call, the PCS phone will emit a tone to alert the user that a nurse call is placed on the system.