

CUSTOMIZABLE SOLUTION

Provider[®] 680
Nurse Call System

Provider 680 Nurse Call System with integrations is a complete patient/resident to staff alerting and communications solution for a nursing unit, a single building facility, or a multi-building campus.

The integrated Provider 680 solution ensures direct alerting and communication for a safer environment, which drives both patient and resident satisfaction.

Nurse Station Options

Touchscreen Nurse Console

- Graphic display of multiple calls
- Embedded technology designed for continuous operation
- Wall or desk-mount option
- Menu-driven touchscreen navigation
- One-touch operation to alert staff
- Communicate with associated stations
- Tone-Visual option



Satellite Nurse Console

- Full functionality in a space saving design
- Wall or desk-mount
- Multiple call-in display



Annunciator Panel

- Dedicated call indications for each room
- Full English call display
- Wall or desk-mount



Staff Area Options

Duty Stations

- Tone/visual indication of outstanding calls
- Tone mute (regenerates with next call)
- Intercom option (to Nurse Consoles)



Patient Stations

- Built-in LED night light
- Normal, Priority, Personal Attention, and Cord Out call priorities
- Ultra-reliable tilt release connector
- Compatible with most bed siderail controls
- Entertainment and light control interface
- Tone-visual station options
- No dummy plug required



Pillow Speakers and Call Cords

- Completely sealed – easily cleaned
- UL listed for use in oxygen-enriched atmosphere
- Large easy-to-press buttons



Bath and Special Function Stations

- Built-in LED night light
- One-touch/pull operation
- Auxiliary alarm input (vent alarm, bed exit, IV alarm)
- Audio to the bathroom



OPTIONS & INTEGRATIONS

SIP WIRELESS PHONES



CORDLESS PHONES



TEXT MESSAGING & POCKET PAGING




AUTOMATED OVERHEAD PAGING



Provider® 680 CUSTOMIZED FOR YOUR WORKFLOWS



REAL-TIME LOCATING



PC CONSOLE ACTIVITY DISPLAY



CALL ACTIVITY REPORTING

Provider® 790 Nurse Call System

CROSS PLATFORM INTEGRATION

Learn more about the Provider 680 Nurse Call System; contact Jeron at **800-621-1903** or www.jeron.com